


Pe Rabbit Anti-Human Enpp3/Cd203C Mab - O14638

Item Number: CM0011885



Introduction

PE-conjugated rabbit monoclonal antibody targeting human ENPP3/CD203c (O14638) for flow cytometry. Ideal for immunophenotyping basophils and mast cells in allergic inflammation research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human ENPP3/CD203c mAb
Remarks/Alias	B10; NPP3; PDNP3; CD203c; PD-IBETA
Species	Human
GeneID (human)	5169
GeneID	5169
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O14638
Protein Weight	100kDa
Shipping	Ice bag