


Pe Rabbit Anti-Human Fgfr3 Monoclonal Antibody For Flow Cytometry - P22607

Item Number: CM0009839



Introduction

PE-conjugated rabbit monoclonal antibody targeting human FGFR3 (CD333) for flow cytometry. Validated for human FC; recombinant immunogen amino acids 23-375. Ideal for studying FGFR3 signaling in cancer and skeletal development.

[Learn More](#)

Parameter	Specification
Product Name	PE Rabbit anti-Human FGFR3 mAb
Remarks/Alias	ACH; CEK2; JTK4; CD333; HSFGR3EX
Species	Human
GeneID (Human)	2261
GeneID	2261
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-375 of human FGFR3 (NP_000133.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P22607
Protein Weight	88kDa
Shipping	Ice bag