


Pe Rabbit Anti-Human Gprc5D Monoclonal Antibody For Flow Cytometry - Q9Nzd1

Item Number: CM0011152



Introduction

PE-conjugated rabbit monoclonal antibody against human GPRC5D for flow cytometry. GPRC5D is a G protein-coupled receptor involved in hard keratin expression and hair/nail development. Suitable for human samples.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human GPRC5D mAb
Remarks/Alias	GPRC5D
Species	Human
Gene ID	55507
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human GPRC5D (NP_061124.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NZD1
Protein Weight	34kDa/36kDa/38kDa
Shipping	Ice bag