

Pe Rabbit Anti-Human Hif-1 Alpha Monoclonal Antibody For Fc (Intra) - Q16665

Item Number: CM0004182



Introduction

PE Rabbit anti-Human HIF-1 alpha Monoclonal Antibody for intracellular flow cytometry. Suitable for hypoxia response, cancer, and metabolic adaptation research. Reacts with human. - Q16665

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human HIF-1 alpha mAb
Remarks/Alias	HIF1, MOP1, PASD8, HIF-1A, bHLHe78, HIF-1alpha, HIF1-ALPHA, HIF-1-alpha
Species	Human
GeneID (Human)	3091
GeneID	3091
Immunogen	Synthetic peptide corresponding to amino acids 401-500 of HIF1 alpha (NP_001521.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	Q16665
Protein Weight	93kDa
Shipping	Ice bag