

Pe Rabbit Anti-Human Hmox1 Mab - P09601

Item Number: CM0009840



Introduction

PE-conjugated rabbit monoclonal antibody against human Heme Oxygenase 1 (HMOX1/HO-1) for intracellular flow cytometry. Recognizes human HMOX1; immunogen corresponds to amino acids 1-265. Suitable for apoptosis and oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human HMOX1 mAb
Remarks/Alias	HO-1; HSP32; HMOX1D; bK286B10
Species	Human
GeneID (Human)	3162
GeneID	3162
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-265 of human Heme Oxygenase 1 (HO-1/HMOX1) (NP_002124.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P09601
Protein Weight	33kDa
Shipping	Ice bag