


Pe Rabbit Anti-Human Icam-1/Cd54 Monoclonal Antibody For Fc - P05362

Item Number: CM0010614



Introduction

High-performance PE-conjugated rabbit monoclonal antibody against human ICAM-1/CD54 (SWISS: P05362) validated for flow cytometry (FC) applications, suitable for immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human ICAM-1/CD54 mAb
Remarks/Alias	BB2; CD54; P3.58
Species	Human
GeneID (Human)	3383
GeneID	3383
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-480 of human ICAM-1/CD54 (NP_000192.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P05362
Protein Weight	58kDa
Shipping	Ice bag