

Pe Rabbit Anti-Human Ifn-Gamma Monoclonal Antibody For Fc (Intra) - P01579

Item Number: CM0012546



Introduction

PE-conjugated rabbit monoclonal antibody targeting human Interferon gamma (IFN-gamma) for intracellular flow cytometry. Ideal for detecting cytokine expression in activated T cells for immune research and diagnostic development.

[Learn More](#)

Parameter	Specification
Product Name	PE Rabbit anti-Human IFN-gamma mAb
Remarks/Alias	IFG; IFI; IMD69
Species	Human
Gene ID	3458 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-166 of human IFN-gamma (NP_000610.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P01579
Protein Weight	19kDa
Shipping	Ice bag