


Pe Rabbit Anti-Human IgA1 Monoclonal Antibody For Flow Cytometry - P01876

Item Number: CM0011612



Introduction

PE-conjugated rabbit monoclonal antibody targeting human IgA1 constant region for flow cytometry. Recombinant immunogen (aa 99-353) ensures specificity for the IgA1 subclass; suitable for human sample analysis.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human IgA1 mAb
Remarks/Alias	IgA1
Species	Human
GeneID (Human)	3493
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 99-353 of human IgA1.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01876
Protein Weight	38kDa/42KDa
Shipping	Ice bag