


Pe Rabbit Anti-Human IgG (Fc) Monoclonal Antibody For Flow Cytometry - P01857

Item Number: CM0012685



Introduction

PE-conjugated rabbit monoclonal antibody targeting human IgG Fc region (IGHG1) for flow cytometry. Recognizes human IgG1 constant region. Recombinant Fc fragment immunogen. Suitable for human samples.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human IgG (Fc) mAb
Remarks/Alias	IGHG1
Species	Human
GeneID (Human)	3500
GeneID	3500
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-330 of human IgG (Fc) (P01857)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01857
Protein Weight	36kDa
Shipping	Ice bag