


Pe Rabbit Anti-Human IgG (Fc) Monoclonal Antibody For Flow Cytometry - P01857

Item Number: CM0011248



Introduction

PE-conjugated rabbit monoclonal antibody targeting the Fc region of human IgG heavy constant gamma 1 (IGHG1). Reacts with human IgG1, suitable for flow cytometry (FC) applications. SWISS: P01857.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human IgG (Fc) mAb
Remarks/Alias	IGHG1
Species	Human
GeneID (Human)	3500
GeneID	3500
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01857
Protein Weight	36 kDa
Shipping	Ice bag