

Pe Rabbit Anti-Human IL-10Ra/Cd210 Monoclonal Antibody For Flow Cytometry - Q13651

Item Number: CM0012155



Introduction

PE-conjugated rabbit monoclonal antibody against human IL-10RA/CD210 (IL-10 receptor alpha), validated for flow cytometry. Ideal for studying IL-10-mediated anti-inflammatory signaling in hematopoietic cells.

[Learn More](#)

Parameter	Description
Product Name	PE Rabbit anti-Human IL-10RA/CD210 mAb
Remarks/Alias	CD210; IL10R; CD210a; CDW210A; HIL-10R; IL-10R1
Species	Human
Gene ID (Human)	3587
Gene ID	3587
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-235 of human IL-10RA/CD210 (NP_001549.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Human
SWISS	Q13651
Protein Weight	63 kDa
Shipping	Ice bag