

Pe Rabbit Anti-Human IL-17Ra/Cd217 Mab For Fc - Q96F46

Item Number: CM0009832



Introduction

PE-conjugated rabbit monoclonal antibody targeting human IL-17RA/CD217, validated for flow cytometry. Ideal for studying IL-17 signaling in innate and adaptive immunity, inflammation, and host-pathogen interactions.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human IL-17RA/CD217 mAb
Remarks/Alias	CD217; IL17R; IMD51; CANDF5; CDw217; IL-17RA; hIL-17R
Species	Human
GeneID (Human)	23765
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-320 of human IL-17RA/CD217 (NP_055154.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q96F46
Protein Weight	96kDa
Shipping	Ice bag