


Pe Rabbit Anti-Human IL-2 Monoclonal Antibody For Flow Cytometry (Intracellular) - P60568

Item Number: CM0008959



Introduction

PE-conjugated rabbit monoclonal antibody against human Interleukin-2 (IL-2), validated for intracellular flow cytometry. Detects human IL-2, a key cytokine regulating T-cell growth and immune tolerance.

[Learn More](#)

Parameter	Details
Product Name	PE Rabbit anti-Human IL-2 mAb
Remarks/Alias	IL-2; TCGF; lymphokine
Species	Human
GeneID (Human)	3558
GeneID	3558
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-153 of human IL-2 (NP_000577.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P60568
Protein Weight	18kDa
Shipping	Ice bag