


## Pe Rabbit Anti-Human IL-4 Mab - P05112

Item Number: CM0007896



### Introduction

PE-conjugated rabbit monoclonal antibody targeting human interleukin-4 (IL-4), validated for intracellular flow cytometry. Suitable for immune function research, including IgE regulation and Th2 response studies.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PE Rabbit anti-Human IL-4 mAb                                                                                                         |
| Remarks/Alias    | BSF1; IL-4; BCGF1; BSF-1; BCGF-1                                                                                                      |
| Species          | Human                                                                                                                                 |
| Gene ID (Human)  | 3565                                                                                                                                  |
| Gene ID          | 3565                                                                                                                                  |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-153 of human IL-4 (NP_000580.1). |
| Source           | Rabbit                                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                                 |
| Application      | FC (intra)                                                                                                                            |
| Cross Reactivity | Human                                                                                                                                 |
| SWISS            | P05112                                                                                                                                |
| Protein Weight   | 17 kDa                                                                                                                                |
| Shipping         | Ice bag                                                                                                                               |