


Pe Rabbit Anti-Human IL-5 Monoclonal Antibody For Flow Cytometry - P05113

Item Number: CM0005683



Introduction

PE-conjugated rabbit anti-human IL-5 monoclonal antibody designed for intracellular flow cytometry (FC). Suitable for human samples. Immunogen: amino acids 20-134 of human IL-5. Ideal for immunology and cytokine research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human IL-5 mAb
Remarks/Alias	EDF, IL-5, TRF, interleukin-5
Species	Human
Gene ID	3567
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 20-134 of human IL-5 (NP_000870.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P05113
Protein Weight	15kDa
Shipping	Ice bag