

Pe Rabbit Anti-Human IL18R1/IL-18 R Alpha Monoclonal Antibody For Fc - Q13478

Item Number: CM0005048



Introduction

PE-conjugated rabbit monoclonal antibody against human IL18R1 (IL-18R alpha, CD218a), validated for flow cytometry. Target molecular weight 62kDa; involved in IL-18-mediated Th1 immune responses.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human IL18R1/IL-18 R alpha mAb
Remarks/Alias	CD218a; IL18RA; IL1RRP; CDw218a; IL-1Rrp; IL-18Ralpha; IL18Ralpha2; IL-18R-alpha
Species	Human
GeneID (Human)	8809
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-329 of human IL18R1/IL-18 R alpha (NP_003846.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13478
Protein Weight	62kDa
Shipping	Ice bag