

Pe Rabbit Anti-Human IL9 Mab - P15248

Item Number: CM0008835



Introduction

PE-conjugated rabbit monoclonal antibody targeting human Interleukin-9 (IL-9, P40/HP40). Validated for intracellular flow cytometry (FC). Ideal for immunology research, detecting cytokine expression in human T-cells and mast cells.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human IL9 mAb
Remarks/Alias	P40; HP40; IL-9
Species	Human
GeneID (human)	3578
GeneID	3578
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 19-144 of human IL9 (NP_000581.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P15248
Protein Weight	16kDa
Shipping	Ice bag