


Pe Rabbit Anti-Human Ilt5/Cd85A Monoclonal Antibody For Fc - O75022

Item Number: CM0009949



Introduction

PE-conjugated rabbit monoclonal antibody targeting human ILT5/CD85a (LILRB3) for flow cytometry. Validated for human cross-reactivity; detects ~69/71 kDa protein.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human ILT5/CD85a mAb
Remarks/Alias	HL9; ILT5; LIR3; PIRB; CD85A; ILT-5; LIR-3; PIR-B; LILRA6
Species	Human
Gene ID (Human)	11025
Gene ID	11025
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O75022
Protein Weight	69kDa/71kDa
Shipping	Ice bag