

Pe Rabbit Anti-Human Integrin Alpha 3/Cd49C Monoclonal Antibody For Flow Cytometry - P26006

Item Number: CM0005099

Introduction

PE-conjugated rabbit monoclonal antibody targeting human Integrin alpha-3/CD49c for flow cytometry applications. Useful in studying cell-extracellular matrix interactions, invadopodia formation, and microbial infection mechanisms.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human Integrin alpha 3/CD49c mAb
Remarks/Alias	JEB7; VL3A; CD49C; FRP-2; GAPB3; ILNEB; MSK18; VCA-2; VLA3a; GAP-B3
Species	Human
Gene ID (Human)	3675
Gene ID	3675
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 466-620 of human Integrin alpha 3/CD49c (NP_002195.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P26006
Protein Weight	117kDa
Shipping	Ice bag