


Pe Rabbit Anti-Human Integrin B5/Itgb5 Monoclonal Antibody For Fc - P18084

Item Number: CM0011440



Introduction

PE-conjugated rabbit monoclonal antibody targeting human Integrin beta-5 (ITGB5) for flow cytometry. Recognizes RGD-dependent fibronectin receptor and adenovirus receptor; cell membrane localization.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human Integrin β 5/ITGB5 mAb
Remarks/Alias	ITGB5; integrin beta-5
Species	Human
GeneID (Human)	3693
GeneID	3693
Immunogen	Recombinant protein corresponding to a sequence within amino acids 24-490 of human Integrin β 5/ITGB5 (NP_002204.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P18084
Protein Weight	88kDa
Shipping	Ice bag