

Pe Rabbit Anti-Human Kir2DL4/CD158D Monoclonal Antibody For Flow Cytometry - Q99706

Item Number: CM0004823



Introduction

PE-conjugated rabbit monoclonal antibody targeting human KIR2DL4 (CD158d), validated for flow cytometry. Recognizes the killer cell immunoglobulin-like receptor involved in NK cell effector functions and maternal-fetal immune tolerance. Cross-reacts with human. - Q99706

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human KIR2DL4/CD158d mAb
Remarks/Alias	G9P; CD158D; KIR103; KIR-2DL4; KIR103AS; KIR-103AS
Species	Human
Gene ID	3805
Immunogen	Recombinant fusion protein containing amino acids 24-242 of human KIR2DL4/CD158d (NP_002246.5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q99706
Protein Weight	41kDa
Shipping	Ice bag