


Pe Rabbit Anti-Human Kir3DL1/Cd158E Monoclonal Antibody For Fc - P43629

Item Number: CM0011651



Introduction

PE-conjugated rabbit monoclonal antibody targeting human KIR3DL1/CD158e (SWISS: P43629). Validated for flow cytometry (FC), detects inhibitory NK cell receptor. Species reactivity: human. Ideal for immunophenotyping.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human KIR3DL1/CD158e mAb
Remarks/Alias	KIR; NKB1; NKAT3; NKB1B; NKAT-3; CD158E1; KIR2DL5B; KIR3DL1/S1
Species	Human
GeneID	3811
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-340 of human KIR3DL1/CD158e (NP_037421.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P43629
Protein Weight	49kDa
Shipping	Ice bag