

Pe Rabbit Anti-Human Lambda Light Chain Monoclonal Antibody For Flow Cytometry - P0Cg04

Item Number: CM0008478



Introduction

PE-conjugated rabbit monoclonal antibody targeting human immunoglobulin lambda constant region 1 (IGLC1). Validated for flow cytometry, ideal for B-cell immunophenotyping and immunoglobulin research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human Lambda Light Chain mAb
Species	Human
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-106 of human Lambda Light Chain.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P0CG04
Protein Weight	11kDa
Shipping	Ice bag