


Pe Rabbit Anti-Human Lamp1/Cd107A Monoclonal Antibody For Intracellular Flow Cytometry - P11279

Item Number: CM0010330



Introduction

PE-conjugated rabbit monoclonal antibody targeting human LAMP1 (CD107a, LAMPA, LGP120) validated for intracellular flow cytometry. Detects lysosomal membrane glycoprotein involved in lysosome biogenesis, autophagy, and NK cell cytotoxicity.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human LAMP1/CD107a mAb
Remarks/Alias	LAMPA; CD107a; LGP120
Species	Human
GeneID (Human)	3916
GeneID	3916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-382 of human LAMP1/CD107a (NP_005552.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P11279
Protein Weight	38kDa/45kDa
Shipping	Ice bag