

Pe Rabbit Anti-Human Lgr5 (Gpr49) Monoclonal Antibody For Flow Cytometry - 075473

Item Number: CM0008654



Introduction

PE-conjugated rabbit monoclonal antibody targeting human LGR5 (GPR49) for flow cytometry. High specificity to C-terminal peptide, ideal for stem cell and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human LGR5 (GPR49) mAb
Remarks/Alias	FEX; HG38; GPR49; GPR67; GRP49
Species	Human
GeneID (human)	8549
GeneID	8549
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 480-580 of human LGR5 (NP_003658.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	075473
Protein Weight	100kDa
Shipping	Ice bag