

Pe Rabbit Anti-Human Lifr/Cd118 Monoclonal Antibody For Fc - P42702

Item Number: CM0004276



Introduction

PE-conjugated rabbit monoclonal antibody targeting human LIFR/CD118 for flow cytometry. Suitable for studying leukemia inhibitory factor receptor signaling and cytokine receptor interactions in immunology and stem cell research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human LIFR/CD118 mAb
Remarks/Alias	SWS; SJS2; STWS; CD118; LIF-R
Species	Human
GeneID (Human)	3977
GeneID	3977
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 45-833 of human LIFR/CD118 (NP_002301.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P42702
Protein Weight	124kDa
Shipping	Ice bag