


Pe Rabbit Anti-Human Mertk Monoclonal Antibody For Flow Cytometry - Q12866

Item Number: CM0011392



Introduction

PE-conjugated rabbit monoclonal antibody targeting human MERTK, validated for flow cytometry. Specific for human. Ideal for studying efferocytosis, cancer signaling, and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human MERTK mAb
Remarks/Alias	MERTK; MER; RP38; Tyro12; c-Eyk; c-mer; tyrosine-protein kinase Mer
Species	Human
GeneID (human)	10461
GeneID	10461
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 26-499 of human MERTK (NP_006334.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q12866
Protein Weight	110kDa
Shipping	Ice bag