

Pe Rabbit Anti-Human Mglur2 Monoclonal Antibody For Fc - Q14416

Item Number: CM0012245



Introduction

PE-conjugated rabbit monoclonal antibody against human mGluR2 (metabotropic glutamate receptor 2), validated for flow cytometry (FC). Ideal for neuroscience, synaptic transmission, and viral entry mechanism studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human mGluR2 mAb
Remarks/Alias	GLUR2; mGlu2; GPRC1B; MGLUR2
Species	Human
GeneID (Human)	2912
GeneID	2912
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-498 of human mGluR2 (NP_000830.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q14416
Protein Weight	96kDa
Shipping	Ice bag