

Pe Rabbit Anti-Human Mica/Micb Monoclonal Antibody For Fc - Q29983 / Q29980

Item Number: CM0005309



Introduction

PE-conjugated rabbit monoclonal antibody against human MICA/MICB for flow cytometry. Reacts with stress-induced NKG2D ligands important for immune surveillance and tumor immunology.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human MICA/MICB mAb
Remarks/Alias	MIC-A; PERB11.1; PERB11.2; MICB
Species	Human
GeneID (human)	100507436/4277
GeneID	100507436/4277
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-298 of human MICB (NP_005922.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q29983/Q29980
Protein Weight	4kDa/31kDa/37kDa/43kDa
Shipping	Ice bag