

Pe Rabbit Anti-Human Mrc2 Monoclonal Antibody For Flow Cytometry - Q9Ubg0

Item Number: CM0009245



Introduction

PE-conjugated rabbit monoclonal antibody targeting human MRC2 (CD280, ENDO180) for flow cytometry. Ideal for studying collagen remodeling, endocytosis, and cancer progression in human cell samples.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human MRC2 mAb
Remarks/Alias	CD280; UPARAP; CLEC13E; ENDO180
Species	Human
Gene ID (Human)	9902
Gene ID	9902
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 35-513 of human MRC2 (NP_006030.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9UBG0
Protein Weight	167kDa
Shipping	Ice bag