

## Pe Rabbit Anti-Human Nectin-1/Cd111 Monoclonal Antibody For Flow Cytometry - Q15223

Item Number: CM0009417



### Introduction

PE-conjugated rabbit monoclonal antibody against human Nectin-1/CD111 (Q15223). Suitable for flow cytometry. Useful for research on cell adhesion, neuronal development, and herpesvirus entry.

[Learn More](#)

| Parameter        | Value                                                                                                                      |
|------------------|----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PE Rabbit anti-Human Nectin-1/CD111 mAb                                                                                    |
| Remarks/Alias    | ED4; PRR; HlgR; HV1S; HVEC; OFC7; PRR1; PVRR; CD111; PVRL1; PVRR1; SK-12; CLPED1; nectin-1                                 |
| Species          | Human                                                                                                                      |
| Gene ID (Human)  | 5818                                                                                                                       |
| Gene ID          | 5818                                                                                                                       |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-334 of human Nectin-1/CD111 (NP_002846.3). |
| Source           | Rabbit                                                                                                                     |
| Category         | Monoclonal Antibodies                                                                                                      |
| Application      | FC                                                                                                                         |
| Cross Reactivity | Human                                                                                                                      |
| SWISS            | Q15223                                                                                                                     |
| Protein Weight   | 57kDa/51kDa/39kDa                                                                                                          |
| Shipping         | Ice bag                                                                                                                    |