


Pe Rabbit Anti-Human Neuropilin-2/Nrp2 Monoclonal Antibody For Fc - O60462

Item Number: CM0009643



Introduction

PE-conjugated rabbit monoclonal antibody targeting human Neuropilin-2 (NRP2) for flow cytometry. High-affinity receptor for semaphorins, VEGF, and HCMV; validated for FC application.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human Neuropilin-2/NRP2 mAb
Remarks/Alias	NP2; NPN2; PRO2714; VEGF165R2
Species	Human
GeneID	8828
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-855 of human Neuropilin-2/NRP2 (NP_957718.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O60462
Protein Weight	105kDa
Shipping	Ice bag