

Pe Rabbit Anti-Human Ng2 Mab For Flow Cytometry - Q6Uvk1

Item Number: CM0008028



Introduction

PE-conjugated rabbit monoclonal antibody targeting human NG2 (CSPG4) for flow cytometry. Recognizes chondroitin sulfate proteoglycan 4, a cell surface receptor involved in cell migration and angiogenesis. Suitable for melanoma and tumor microenvironment studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human NG2 mAb
Remarks/Alias	NG2; MCSP; MCSPG; MSK16; CSPG4A; HMW-MAA; MEL-CSPG
Species	Human
GeneID (Huma)	1464
GeneID	1464
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1583-2224 of human NG2 (NP_001888.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q6UVK1
Protein Weight	251kDa
Shipping	Ice bag