

Pe Rabbit Anti-Human Nt5E/Cd73 Monoclonal Antibody For Flow Cytometry - P21589

Item Number: CM0007730



Introduction

PE-conjugated rabbit monoclonal antibody targeting human NT5E/CD73, validated for flow cytometry. Recognizes ecto-5'-nucleotidase, a GPI-anchored enzyme involved in nucleotide metabolism and immunosuppression. Ideal for immunology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human NT5E/CD73 mAb
Remarks/Alias	NT; eN; NT5; NTE; eNT; CD73; E5NT; CALJA
Species	Human
Gene ID (Human)	4907
Gene ID	4907
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-547 of human NT5E/CD73 (NP_002517.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P21589
Protein Weight	63kDa
Shipping	Ice bag