

## Pe Rabbit Anti-Human Pd-L1/Cd274 Monoclonal Antibody For Fc - Q9Nzq7

Item Number: CM0010810



### Introduction

PE-conjugated rabbit monoclonal antibody against human PD-L1/CD274 (B7-H1, PDL1). Suitable for flow cytometry (FC). Key tool for studying immune tolerance, checkpoint blockade, and cancer immunotherapy. Recognizes a 33kDa protein.

[Learn More](#)

| Parameter        | Value                                              |
|------------------|----------------------------------------------------|
| Product Name     | PE Rabbit anti-Human PD-L1/CD274 mAb               |
| Remarks/Alias    | B7-H; B7H1; PDL1; PD-L1; hPD-L1; PDCD1L1; PDCD1LG1 |
| Species          | Human                                              |
| GeneID (Human)   | 29126                                              |
| GeneID           | 29126                                              |
| Immunogen        | Recombinant protein                                |
| Source           | Rabbit                                             |
| Category         | Monoclonal Antibodies                              |
| Application      | FC                                                 |
| Cross Reactivity | Human                                              |
| SWISS            | Q9NZQ7                                             |
| Protein Weight   | 33kDa                                              |
| Shipping         | Ice bag                                            |