


Pe Rabbit Anti-Human Phospho-Histone H2Ax-S139 Monoclonal Antibody For Flow Cytometry - P16104

Item Number: CM0001699



Introduction

PE-conjugated rabbit monoclonal antibody targeting phospho-S139 of human Histone H2AX, validated for intracellular flow cytometry. Suitable for DNA damage and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human Phospho-Histone H2AX-S139 mAb
Remarks/Alias	H2A.X; H2A/X; H2AFX
Species	Human
GeneID (Human)	3014
GeneID	3014
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P16104
Protein Weight	15kDa
Shipping	Ice bag