

Pe Rabbit Anti-Human Plzf Monoclonal Antibody For Fc (Intra) - Q05516

Item Number: CM0009959



Introduction

PE-conjugated rabbit monoclonal antibody targeting human PLZF (ZBTB16) for intracellular flow cytometry. Validated in human samples, 74kDa, suitable for hematopoietic and transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human PLZF mAb
Remarks/Alias	PLZF; ZNF145
Species	Human
GeneID (Human)	7704
GeneID	7704
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 6-137 of human PLZF (NP_005997.2))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	Q05516
Protein Weight	74kDa
Shipping	Ice bag