


Pe Rabbit Anti-Human Podoplanin Monoclonal Antibody For Flow Cytometry - Q86YI7

Item Number: CM0010123



Introduction

PE-conjugated rabbit monoclonal antibody targeting human podoplanin (PDPN), validated for flow cytometry. Useful for investigating cell migration, adhesion, and epithelial-mesenchymal transition. SWISS: Q86YL7.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human Podoplanin mAb
Remarks/Alias	T1A; GP36; GP40; Gp38; OTS8; T1A2; T11A; D2-40; T1A-2; AGGRUS; HT1A-1; PA2.26
Species	Human
GeneID	10630
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-123 of Human Podoplanin (NP_006465.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q86YL7
Protein Weight	17kDa
Shipping	Ice bag