


Pe Rabbit Anti-Human Pth1R Monoclonal Antibody For Fc - Q03431

Item Number: CM0012153



Introduction

PE-conjugated rabbit monoclonal antibody targeting human PTH1R (parathyroid hormone 1 receptor) for flow cytometry. Recognizes extracellular domain (aa 27-188). Reacts with human PTH1R. ~66 kDa. Research use.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human PTH1R mAb
Remarks/Alias	PFE, EKNS, PTHR, PTHR1
Species	Human
GeneID	5745
Immunogen	Recombinant protein, recombinant fusion protein containing a sequence corresponding to amino acids 27-188 of human PTH1R (NP_000307.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q03431
Protein Weight	66kDa
Shipping	Ice bag