


Pe Rabbit Anti-Human Rankl/Cd254 Monoclonal Antibody For Flow Cytometry (Intracellular) - O14788

Item Number: CM0012196



Introduction

PE-conjugated rabbit monoclonal antibody targeting human RANKL/CD254 (TNFSF11) for intracellular flow cytometry. Ideal for studying osteoclast differentiation, bone remodeling, and dendritic cell/T-cell interactions.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human RANKL/CD254 mAb
Remarks/Alias	ODF; OPGL; sOdf; CD254; OPTB2; RANKL; TNLG6B; TRANCE; hRANKL2
Species	Human
Gene ID	8600
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 136-317 of human RANKL/CD254 (NP_003692.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	O14788
Protein Weight	27kDa/30kDa/35kDa
Shipping	Ice bag