


Pe Rabbit Anti-Human Siglec-7/Cd328 Monoclonal Antibody For Fc - Q9Y286

Item Number: CM0006162



Introduction

PE-conjugated rabbit monoclonal antibody against human Siglec-7 (CD328/AIRM1/p75). Specifically designed for flow cytometry (FC). Recombinant immunogen covering aa 18-144. Ideal for natural killer cell and immune research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human Siglec-7/CD328 mAb
Remarks/Alias	p75; QA79; AIRM1; CD328; AIRM-1; CDw328; D-siglec; SIGLEC-7; SIGLECP2; SIGLEC19P; p75/AIRM1
Species	Human
GeneID (human)	27036
GeneID	27036
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 18-144 of human Siglec-7/CD328 (NP_055200.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9Y286
Protein Weight	51kDa
Shipping	Ice bag