

Pe Rabbit Anti-Human Siglec-8 Monoclonal Antibody For Fc - Q9Nyz4

Item Number: CM0006216



Introduction

PE-conjugated rabbit monoclonal antibody targeting human Siglec-8, validated for flow cytometry (FC). Detects multiple isoforms (44/46/54 kDa) and is ideal for studying eosinophil, basophil, and mast cell biology in human samples.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human Siglec-8 mAb
Remarks/Alias	SAF2; SIGLEC-8; SIGLEC8L
Species	Human
GeneID (Human)	27181
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-363 of human Siglec-8 (NP_055257.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NYZ4
Protein Weight	44kDa/46kDa/54kDa
Shipping	Ice bag