


Pe Rabbit Anti-Human Stat1 Monoclonal Antibody For Flow Cytometry - P42224

Item Number: CM0012164



Introduction

PE-conjugated rabbit monoclonal antibody targeting human STAT1. Validated for intracellular flow cytometry. Ideal for interferon signaling, antiviral defense, and STAT1 transcription factor studies. Research-use IVD raw material.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human STAT1 mAb
Remarks/Alias	CANDF7; IMD31A; IMD31B; IMD31C; ISGF-3; STAT91
Species	Human
GeneID (Human)	6772
GeneID	6772
Immunogen	Synthetic peptide corresponding to amino acids 650-750 of human STAT1 (NP_009330.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P42224
Protein Weight	83 kDa/87 kDa
Shipping	Ice bag