

Pe Rabbit Anti-Human Stat3 Monoclonal Antibody For Intracellular Flow Cytometry - P40763

Item Number: CM0010525



Introduction

PE-conjugated rabbit anti-human STAT3 monoclonal antibody for intracellular flow cytometry (FC). Targets STAT3/APRF, an 88 kDa transcription activator in IL-6 and JAK/STAT pathways; cross-reacts with human. Research use only.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human STAT3 mAb
Remarks/Alias	APRF; HIES; ADMIO; ADMIO1
Species	Human
GeneID (Human)	6774
GeneID	6774
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 611-710 of human STAT3 (NP_003141.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P40763
Protein Weight	88 kDa
Shipping	Ice bag