


Pe Rabbit Anti-Human Tie-2/Cd202B Monoclonal Antibody For Flow Cytometry - Q02763

Item Number: CM0008773



Introduction

PE-conjugated rabbit anti-human Tie-2/CD202b monoclonal antibody specifically validated for flow cytometry (FC) detection of Tie-2 receptor (UniProt Q02763) in human samples. Ideal for angiogenesis and vascular biology research applications.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human Tie-2/CD202b mAb
Remarks/Alias	TIE2; VMCM; GLC3E; TIE-2; VMCM1; CD202B
Species	Human
Gene ID	7010
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-745 of human Tie-2/CD202b (NP_000450.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q02763
Protein Weight	126kDa
Shipping	Ice bag