


Pe Rabbit Anti-Human Tigit Monoclonal Antibody For Flow Cytometry - Q495A1

Item Number: CM0008222



Introduction

PE-conjugated rabbit monoclonal antibody against human TIGIT (VSIG9, VSTM3), validated for flow cytometry. Targets an inhibitory immunoreceptor on T cells and NK cells. UniProt Q495A1.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human TIGIT mAb
Remarks/Alias	VSIG9; VSTM3; WUCAM
Species	Human
GeneID (Human)	201633
GeneID	201633
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q495A1
Protein Weight	18kDa/26kDa
Shipping	Ice bag