


Pe Rabbit Anti-Human Tim-3/Havcr2 Monoclonal Antibody For Fc - Q8Tdq0

Item Number: CM0008423



Introduction

Rabbit monoclonal PE-conjugated antibody against human TIM-3/HAVCR2, validated for flow cytometry. Recombinant immunogen (aa 22-200). Ideal for immune checkpoint and T-cell exhaustion research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human TIM-3/HAVCR2 mAb
Remarks/Alias	TIM3; CD366; KIM-3; SPTCL; TIMD3; Tim-3; TIMD-3; HAVcr-2
Species	Human
Gene ID (Human)	84868
Gene ID	84868
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-200 of human TIM-3/HAVCR2 (NP_116171.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8TDQ0
Protein Weight	33kDa
Shipping	Ice bag