


Pe Rabbit Anti-Human Tlr3 Monoclonal Antibody For Flow Cytometry (Intracellular) - O15455

Item Number: CM0008254



Introduction

PE-conjugated rabbit monoclonal antibody targeting human TLR3 (CD283), validated for intracellular flow cytometry. Ideal for innate immunity and viral infection research applications.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human TLR3 mAb
Remarks/Alias	CD283; IIAE2; IMD83
Species	Human
Gene ID (Human)	7098
Gene ID	7098
Immunogen	Recombinant protein corresponding to a sequence within amino acids 24-339 of human TLR3 (NP_003256.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	O15455
Protein Weight	104kDa
Shipping	Ice bag