

Pe Rabbit Anti-Human Tmprss2 Monoclonal Antibody For Flow Cytometry - O15393

Item Number: CM0004793



Introduction

PE-conjugated rabbit anti-human TMPRSS2 monoclonal antibody for flow cytometry. Suitable for studying prostate cancer progression, viral infection mechanisms (SARS-CoV, SARS-CoV-2), and host cell proteolysis.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human TMPRSS2 mAb
Remarks/Alias	PRSS10
Species	Human
Gene ID (Human)	7113
Gene ID	7113
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 106-255 of human TMPRSS2 (NP_005647.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O15393
Protein Weight	54kDa
Shipping	Ice bag