


Pe Rabbit Anti-Human Tnf-A Monoclonal Antibody For Fc (Intra) - P01375

Item Number: CM0007785



Introduction

PE-conjugated rabbit anti-human TNF- α monoclonal antibody for flow cytometry intracellular staining. Investigate tumor necrosis factor alpha in inflammatory responses, immune regulation, and cancer biology. Cross-reacts with human TNF- α (cachectin/TNFSF2).

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human TNF- α mAb
Remarks/Alias	DIF; TNFA; TNFSF2; TNLG1F; TNF-alpha
Species	Human
Gene ID (Human)	7124
Gene ID	7124
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 77-233 of human TNF- α (NP_002344.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P01375
Protein Weight	26kDa
Shipping	Ice bag